



SARS-CoV-2 Antigen Quick Test

Convenient Testing from Saliva

The Test



- ✓ Simple
- ✓ Convenient
- ✓ Fast
- ✓ Sensitive



Performance Summary

Sample material	saliva
Time to result	10 min
Limit of detection	1.44 x 10³ TCID₅₀/ml
Clinical sensitivity	94% (33/35)
Clinical specificity	100% (70/70)

The Kit

[COV04S] PCL COVID19 Ag Gold Saliva

All you need for saliva testing

Available from GFE



 Produced by PCL Inc., Republic of Korea

Contact

Information

sales@gfeblut.de

Order

order@gfeblut.de